

Product Datasheet

F(ab')₂ Goat anti-Rabbit IgG Heavy and Light Chain Cross-Adsorbed Antibody HRP Conjugated (orb1531016)

Description	Rabbit IgG-Heavy and Light Chain Antibody
Species/Host	Goat
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC, WB
Target	Rabbit IgG-Heavy and Light Chain
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.05% Pro-Clean 400
Form/Appearance	F(ab') ₂
Concentration	0.5 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: F(ab') ₂
Clonality	Polyclonal
Dilution Range	WB - 1:5,000 - 1:50,000; IHC - 1:200 - 1:5,000; ICC - 1:200 - 1:5,000; ELISA - 1:10,000 - 1:100,000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)