

## Product Datasheet

### S-Tag/S15 Antibody (orb1532949)

<b>Description</b>	S-Tag/S15 Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Bovine
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IF, WB
<b>Immunogen</b>	Peptide derived from the N-terminal sequence of pET30a(+)
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG
<b>Dilution Range</b>	Immunofluorescence: N/A; Western Blot: 1:2,000-5,000
<b>Expiration Date</b>	12 months from date of receipt.