

Product Datasheet

S-Tag/S15 Antibody (orb1532952)

Description	S-Tag/S15 Antibody
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide derived from the N-terminal sequence of pET30a(+)
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Dilution Range	Western Blot: 1:2,000-5,000
Expiration Date	12 months from date of receipt.