


Product Datasheet

V5 epitope antibody (orb153323)

Description	Goat polyclonal antibody to V5 epitope. V5 peptide corresponds to amino acid residues GKPIPPELLGLDST...
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Purified recombinant peptide produced in E. coli.
Target	V5
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Polyclonal antibody supplied as 200 or 500 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Protein Sequence	MGKPIPPELLGLDSTNPMGKPIPPELLGLDSTNPMGKPIPPELLGLDSTNPMGKPIPPELLGLDSTNPMGKPIPPELLGLDSTNPGRLERPHRD
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:500-1:5,000, IF:1:50-1:2,000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell line
I...