

Product Datasheet

SARS-CoV-2 Nucleoprotein Antibody (orb1533465)

Description	SARS-CoV-2 Nucleoprotein Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	LP-13473, synthetic peptide
Target	SARS-CoV-2 Nucleoprotein
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	ELISA, WB
Expiration Date	12 months from date of receipt.