

## Product Datasheet

### SARS-CoV-2 Nucleoprotein Antibody (orb1533466)

<b>Description</b>	SARS-CoV-2 Nucleoprotein Antibody
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	LP-13474, synthetic peptide
<b>Target</b>	SARS-CoV-2 Nucleoprotein
<b>Preservatives</b>	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	ELISA, WB
<b>Expiration Date</b>	12 months from date of receipt.