



Product Datasheet

PIK3R1 Antibody (orb153366)


Description	Mouse monoclonal antibody to PIK3R1
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC
Immunogen	Purified recombinant fragment of human PIK3R1 (AA: 159-388) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	83.6kDa
Entrez	5295
Dilution Range	IHC: 1/200 - 1/1000; ICC: 1/100; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...



Flow cytometric analysis of NIH3T3 cells...



Immunohistochemical analysis of formalin...