
Product Datasheet

HTR1F / 5-HT1F Receptor Antibody (C-Terminus) (orb1533938)

Description

HTR1F / 5-HT1F Receptor Antibody

Species/Host

Rabbit

Reactivity

Bovine, Canine, Gallus, Guinea pig, Hamster, Human, Mammal, Monkey, Mouse, Rabbit, Rat

Conjugation

Unconjugated

Tested Applications

IHC, IHC-P

Immunogen

Synthetic 17 amino acid peptide from C-Terminus of human 5HT1F Receptor. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Orangutan, Gibbon, Monkey, Mouse, Rat, Hamster, Elephant, Panda, Dog, Bovine, Bat, Rabbit, Guinea pig, Turkey, Chicken, Platypus (100%); Xenopus, Stickleback, Pufferfish (94%); Opossum, Lizard (88%); Horse (82%).

Target

HTR1F / 5-HT1F Receptor

Preservatives

PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide

Concentration

1 mg/ml

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Clonality

Polyclonal

Dilution Range

IHC, IHC-P (2 - 3 µg/ml)

Expiration Date

12 months from date of receipt.