



Product Datasheet

GSC Antibody (orb153469)


Description	Mouse monoclonal antibody to GSC
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of human GSC (AA: 191-257) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	28.2kDa
Entrez	145258
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...



Western blot analysis of human GSC using...



Western blot analysis of HEK293 (Lane 1)...