

Product Datasheet

C14orf1 Antibody (C-Terminus) (orb1534701)

Description	C14orf1 Antibody (C-Terminus)
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mammal, Monkey, Rabbit
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 15 amino acid peptide from C-Terminus of human C14orf1. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Gibbon, Monkey, Marmoset, Elephant, Panda, Bat, Rabbit, Guinea pig (100%); Galago, Mouse, Rat, Horse (93%); Platypus (87%); Dog, Bovine, Pig (80%).
Target	C14orf1
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (10 - 15 µg/ml)
Expiration Date	12 months from date of receipt.