

Product Datasheet

C14orf1 Antibody (Internal) (orb1534703)

Description	C14orf1 Antibody (Internal)
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Hamster, Human, Monkey, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 14 amino acid peptide from internal region of human C14orf1. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Gibbon, Monkey, Galago, Marmoset, Mouse, Rat, Hamster, Elephant, Panda, Dog, Horse, Rabbit, Opossum, Guinea pig (100%); Bovine, Bat, Pig, Platypus, Salmon, Stickleback (93%); Smelt, Medaka, Pufferfish (86%).
Target	C14orf1
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (10 - 15 µg/ml)
Expiration Date	12 months from date of receipt.