

Product Datasheet

SLC12A2/NKCC1 Antibody (orb1534715)

Description	SLC12A2/NKCC1 Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Guinea pig, Human, Monkey, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 15 amino acid peptide from internal region of human SLC12A2 / NKCC1. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Orangutan, Gibbon, Monkey, Marmoset, Mouse, Rat, Ferret, Elephant, Panda, Bovine, Dog, Horse, Pig, Opossum, Guinea pig (100%); Bat (93%).
Target	SLC12A2 / NKCC1
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (5 - 10 µg/ml)
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)