

Product Datasheet

SLC29A1 / ENT1 Antibody (Internal) (orb1534716)

Description	SLC29A1 / ENT1 Antibody (Internal)
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 16 amino acid peptide from internal region of human ENT1. Percent identity with other species by BLAST analysis: Human, Chimpanzee, Gorilla, Orangutan, Gibbon, Monkey, Marmoset (100%); Galago, Horse, Rabbit, Guinea pig (81%).
Target	SLC29A1 / ENT1
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (15 - 20 µg/ml)
Expiration Date	12 months from date of receipt.