

Product Datasheet

SLC31A1 / CTR1 Antibody (Internal) (orb1534791)

Description	SLC31A1 / CTR1 Antibody (Internal)
Species/Host	Rabbit
Reactivity	Canine, Human, Monkey, Rabbit
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 16 amino acid peptide from internal region of human CTR1. Percent identity with other species by BLAST analysis: Human, Gorilla, Orangutan, Gibbon, Monkey, Marmoset, Dog, Elephant, Panda, Rabbit, Opossum, Platypus (100%); Mouse, Rat, Bovine, Hamster, Horse, Pig, Turkey, Chicken (94%); Xenopus (88%); Pufferfish (81%).
Target	SLC31A1 / CTR1
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (15 µg/ml)
Expiration Date	12 months from date of receipt.