

Product Datasheet

USP20 / VDU2 Antibody (Internal) (orb1534810)

Description	USP20 / VDU2 Antibody (Internal)
Species/Host	Rabbit
Reactivity	Canine, Equine, Hamster, Human, Rabbit
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P
Immunogen	Synthetic 18 amino acid peptide from internal region of human USP20. Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Monkey, Hamster, Elephant, Dog, Horse, Rabbit (100%); Orangutan, Marmoset, Mouse, Rat, Panda, Bat, Bovine, Pig (94%); Opossum, Platypus (89%); Turkey, Chicken (83%).
Target	USP20 / VDU2
Preservatives	PBS, 0.1% Sodium Azide. PBS, 0.1% Sodium Azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (10 µg/ml)
Expiration Date	12 months from date of receipt.