
Product Datasheet

**ENPP2 / Autotaxin Antibody (aa698-712)
(orb1535167)**

Description

ENPP2 / Autotaxin Antibody (aa698-712)

Species/Host

Goat

Reactivity

Human, Monkey

Conjugation

Unconjugated

Tested Applications

ELISA, IHC, IHC-P, WB

Immunogen

Peptide with sequence PDHLTSCVRPDVRVS, from the internal region of the protein sequence according to NP_006200.3NP_001035181.1, NP_001124335.1.

Target

ENPP2 / Autotaxin

Preservatives

TBS, pH 7.3, 0.02% Sodium Azide, 0.5% rAlbumin

Concentration

0.5 mg/ml

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

Further information: Peptide ELISA: antibody detection limit dilution 1:128000. Western blot: Approx 100kD band observed in Human Placenta lysates (calculated MW of 105kD according to NP_006200.3). Recommended concentration: 0.3-1 ug/ml.

Clonality

Polyclonal

Dilution Range

IHC, IHC-P (5 µg/ml), Peptide-ELISA (1:128000), WB (0.3 - 1 µg/ml)

Expiration Date

12 months from date of receipt.