

Product Datasheet

GLB1 / Beta-Galactosidase Antibody (clone 10B2) (orb1535248)

Description	GLB1 / Beta-Galactosidase Antibody (clone 10B2)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IHC-P, WB
Immunogen	Full length human recombinant protein of human GLB1(NP_001073279) produced in HEK293T cell.
Target	GLB1 / Beta-Galactosidase
Preservatives	PBS, pH 7.3, 0.02% Sodium Azide, 50% Glycerol, 1% rAlbumin
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IF (1:100), IHC, IHC-P (10 µg/ml), WB (1:200 - 1:2000)
Expiration Date	12 months from date of receipt.