

Product Datasheet

AMY2A / Pancreatic Amylase Antibody (aa484-511) (orb1535268)

Description	AMY2A / Pancreatic Amylase Antibody (aa484-511)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, IHC-P, WB
Target	AMY2A / Pancreatic Amylase
Preservatives	PBS, 0.09% Sodium Azide. PBS, 0.09% Sodium Azide
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IF (1:10 - 1:50), IHC, IHC-P (1:250), WB (1:1000)
Expiration Date	12 months from date of receipt.