

Product Datasheet

PP2CB / PPP2CB Antibody (orb1535301)

Description	PP2CB / PPP2CB Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P, WB
Immunogen	Recombinant human Serine/threonine-protein phosphatase 2A catalytic subunit beta isoform protein.
Target	PP2CB / PPP2CB
Preservatives	PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.. PBS, pH 7.4, 0.03% Proclin 300, 50% glycerol.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Further information: Recommended conditions for Immunohistochemistry: For FFPE tissues, perform steam-based antigen retrieval using Tris/EDTA buffer at pH 9.0.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA, IHC, IHC-P (20 µg/ml), WB
Expiration Date	12 months from date of receipt.