

## Product Datasheet

### VDAC1 / PORIN Antibody (orb1535434)

<b>Description</b>	VDAC1 / PORIN Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, IHC-P, WB
<b>Immunogen</b>	Synthetic peptide conjugated to KLH.
<b>Target</b>	VDAC1 / PORIN
<b>Preservatives</b>	PBS, pH 7.4, 0.02% Sodium Azide. PBS, pH 7.4, 0.02% Sodium Azide
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	ELISA (1:20000), IHC, IHC-P (5 µg/ml), WB (1:250)
<b>Expiration Date</b>	12 months from date of receipt.