

Product Datasheet

COMMD7 Antibody (Internal) (orb1535558)

Description	COMMD7 Antibody (Internal)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the center region of human COMMD7. The exact sequence is proprietary.
Target	COMMD7
Preservatives	0.1 M Tris, pH 7, 0.1 M Glycine, 10% Glycerol, 0.01% Thimerosal. 0.1 M Tris, pH 7, 0.1 M Glycine, 10% Glycerol, 0.01% Thimerosal
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Further information: IHC-P: Suggested antigen retrieval using heat mediated 10 mM Citrate buffer (pH 6.0).
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (5 µg/ml), WB (1:500 - 1:3000)
Expiration Date	12 months from date of receipt.