

## Product Datasheet

# NFAT4 / NFATC3 Antibody (clone 3A12) (orb1535656)

<b>Description</b>	NFAT4 / NFATC3 Antibody (clone 3A12)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC, IHC-P
<b>Immunogen</b>	NFATC3 (NP_775188, 70 aa ~ 150 aa) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
<b>Target</b>	NFAT4 / NFATC3
<b>Preservatives</b>	PBS, pH 7.4. PBS, pH 7.4
<b>Concentration</b>	0.44 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Further information: Immunohistochemistry: Formalin-fixed paraffin-embedded sections. Sandwich ELISA: Recombinant protein.
<b>Isotype</b>	IgG2b,k
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	ELISA, IF (10 µg/ml), IHC, IHC-P (1.5 - 5 µg/ml)
<b>Expiration Date</b>	12 months from date of receipt.