

## Product Datasheet

### OXSR1 / OSR1 Antibody (Internal) (orb1535885)

<b>Description</b>	OXSR1 / OSR1 Antibody (Internal)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Bovine, Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, IHC-P, WB
<b>Immunogen</b>	KLH-conjugated synthetic peptide encompassing a sequence within the central region of human OXSR1.
<b>Target</b>	OXSR1 / OSR1
<b>Preservatives</b>	PBS, pH 7.3, 0.01% Sodium Azide, 30% Glycerol. PBS, pH 7.3, 0.01% Sodium Azide, 30% Glycerol
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	IHC, IHC-P (1:100), WB (1:500 - 1:1000)
<b>Expiration Date</b>	12 months from date of receipt.