

## Product Datasheet

# Tubulin Gamma Antibody (clone TU-32) (orb1536521)

<b>Description</b>	Tubulin Gamma Antibody (clone TU-32)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Animal, Plant, Protozoa
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IHC, IHC-P, WB
<b>Immunogen</b>	Human gamma-tubulin peptide EYHAATRPDYISWGTO, amino acids 434-449.
<b>Target</b>	Tubulin Gamma
<b>Preservatives</b>	PBS, pH 7.4, 15 mM Sodium Azide. PBS, pH 7.4, 15 mM Sodium Azide
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	ICC, IHC, IHC-P (10 µg/ml), WB
<b>Expiration Date</b>	12 months from date of receipt.