

Product Datasheet

PADI4 / PAD4 Antibody (aa1-300) (orb1536575)

Description	PADI4 / PAD4 Antibody (aa1-300)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	Recombinant PADI4 (Met1-Pro300)
Target	PADI4 / PAD4
Preservatives	0.01 M PBS, pH 7.4, 0.05% ProClin 300, 50% Glycerol
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Further information: ELISA 1:100-1:5000; ICC 1:50-500; WB 1:50-400;
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (20 µg/ml), WB
Expiration Date	12 months from date of receipt.