

## Product Datasheet

# GAD 65+67 Antibody (520-600 aa, C-terminal) (orb1536839)

<b>Description</b>	GAD 65+67 Antibody (520-600 aa, C-terminal)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC, IHC-P, WB
<b>Immunogen</b>	Synthesized peptide derived from human GAD-65/67
<b>Target</b>	GAD 65+67
<b>Preservatives</b>	PBS, 0.02% Sodium Azide, 50% Glycerol, 0.5% rAlbumin
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	ELISA (1:5000), IF (1:200 - 1:1000), IHC (1:100 - 1:300), IHC-P (5 µg/ml), WB (1:500 - 1:2000)
<b>Expiration Date</b>	12 months from date of receipt.