

Product Datasheet

PADI4 / PAD4 Antibody (aa2-14) (orb1536867)

Description	PADI4 / PAD4 Antibody (aa2-14)
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Peptide with sequence AQGTLIRVTPEQPC, from the N-Terminus of protein sequence according to NP_036519.2.
Target	PADI4 / PAD4
Preservatives	TBS, pH 7.3, 0.02% Sodium Azide, 0.5% rAlbumin
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC-P (5 µg/ml), Peptide-ELISA (1:64000), WB (0.3 - 1 µg/ml)
Expiration Date	12 months from date of receipt.