

## Product Datasheet

### EBV LMP Antibody (clone CS1-4) (orb1537385)

<b>Description</b>	EBV LMP Antibody (clone CS1-4)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P
<b>Target</b>	EBV LMP
<b>Preservatives</b>	Filter sterilized and diluted in buffer, pH 7.5, rAlbumin, Sodium Azide
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Further information: Immunohistochemistry: Paraffin and frozen.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	IHC-P (1:50)
<b>Expiration Date</b>	12 months from date of receipt.