

Product Datasheet

RPL22/Ribosomal Protein L22 Antibody (orb1538451)

Description	RPL22/Ribosomal Protein L22 Antibody
Species/Host	Goat
Reactivity	Bovine, Canine, Equine, Gallus, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Xenopus
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P
Immunogen	Peptide with sequence C-SKESYELRYFQINQ, from the C-Terminus of protein sequence according to NP_000974.1.
Target	RPL22 / Ribosomal Protein L22
Preservatives	TBS, pH 7.3, 0.02% Sodium Azide, 0.5% rAlbumin
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (3.75 µg/ml), Peptide-ELISA (1:16000)
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)