

---

## Product Datasheet

### **SYNGR3 / Synaptogyrin 3 Antibody (orb1538481)**

**Description**

SYNGR3 / Synaptogyrin 3 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested Applications**

ELISA, IF, IHC, IHC-P, WB

**Immunogen**

A 19 amino acid synthetic peptide near the C-Terminus of human SYNGR3. The immunogen is located within the last 50 amino acids of SYNGR3.

**Target**

SYNGR3 / Synaptogyrin 3

**Preservatives**

PBS containing 0.02% Sodium Azide. PBS containing 0.02% Sodium Azide

**Concentration**

1 mg/ml

**Storage**

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

**Note**

For research use only

**Application notes**

Further information: SYNGR3 antibody can be used for detection of SYNGR3 by Western blot at 1 - 2 ug/mL. Antibody can also be used for immunohistochemistry starting at 2.5 ug/mL. For immunofluorescence start at 20 ug/mL. Antibody validated: Western Blot in mouse samples; Immunohistochemistry in rat samples and Immunofluorescence in rat samples. All other applications and species not yet tested. Western Blot:

**Isotype**

IgG

**Clonality**

Polyclonal

**Dilution Range**

ELISA, IF, IHC, IHC-P (5 µg/ml), WB (1 - 2 µg/ml)

**Expiration Date**

12 months from date of receipt.