
Product Datasheet

NLRP14 / NALP14 Antibody (orb1538544)

Description

NLRP14 / NALP14 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, IHC, IHC-P, WB

Immunogen

A 16 AA synthetic peptide near the N-Terminus of human NALP14. The immunogen is located within aa 200-250 of NALP14.

Target

NLRP14 / NALP14

Preservatives

PBS containing 0.02% Sodium Azide. PBS containing 0.02% Sodium Azide

Concentration

1 mg/ml

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

Further information: NALP14 antibody can be used for detection of NALP14 by Western blot at 1 ug/mL. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested. Western Blot:

Isotype

IgG

Clonality

Polyclonal

Dilution Range

ELISA, IHC, IHC-P (5 µg/ml), WB (1 - 2 µg/ml)

Expiration Date

12 months from date of receipt.