

Product Datasheet

CCL20 / MIP-3-Alpha Antibody (aa33-44) (orb1538673)

Description	CCL20 / MIP-3-Alpha Antibody (aa33-44)
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IHC-P
Immunogen	Peptide with sequence CLGYTDRILHPK, from the internal region of the protein sequence according to NP_004582.1NP_001123518.1.
Target	CCL20 / MIP-3-Alpha
Preservatives	TBS, pH 7.3, 0.02% Sodium Azide, 0.5% rAlbumin
Concentration	0.5 mg/ml
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Further information: Peptide ELISA: antibody detection limit dilution 1:16000. Western blot: Preliminary experiments in Human Peripheral Blood Mononucleocytes lysates gave no specific signal but low background (at antibody concentration up to 2 ug/ml).
Clonality	Polyclonal
Dilution Range	IHC, IHC-P (6.12 µg/ml), Peptide-ELISA (1:16000)
Expiration Date	12 months from date of receipt.