

Product Datasheet

AGO2/EIF2C2 Antibody (orb1538731)


Description	AGO2/EIF2C2 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IHC-P, WB
Immunogen	EIF2C2 (AAH07633, 1 aa ~ 377 aa) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Target	AGO2 / EIF2C2
Preservatives	PBS, pH 7.4. PBS, pH 7.4
Concentration	0.43 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Further information: Immunohistochemistry: Formalin-fixed paraffin-embedded sections. RNAi Knockdown: Antibody validated. Western Blot with the recombinant protein used as the immunogen.
Isotype	IgG1,k
Clonality	Monoclonal
Clone Number	2E12-1C9
Dilution Range	ELISA, IF (10 µg/ml), IHC, IHC-P (3 - 5 µg/ml), WB
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)



Immunofluorescence of monoclonal antibody to EIF2C2 on HeLa cell. [antibody concentration 10 ug/ml]

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44%201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)