



Product Datasheet

Asymmetric DiMethyl-Histone H4-R3 antibody (orb153927)


Description	Rabbit polyclonal antibody to Asymmetric DiMethyl-Histone H4-R3
Reactivity	All, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	A synthetic asymmetric dimethylated peptide around R3 of human histone H4(NP_003539.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	11kDa
Uniprot ID	P62805
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200, IP: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HeLa...



Immunohistochemical staining of paraffin...



IHC-P image of rat kidney tissue using H...