

Product Datasheet

Myoglobin mouse Monoclonal Antibody FITC Conjugated (orb1539557)

Description	Myoglobin mouse Monoclonal Antibody FITC Conjugated
Species/Host	Mouse
Conjugation	FITC
Tested Applications	ICC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	mouse
Uniprot ID	P02144
Dilution Range	ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.