

Product Datasheet

CD16 mouse Monoclonal Antibody HRP Conjugated (orb1539590)

| | |
|----------------------------|---|
| Description | CD16 mouse Monoclonal Antibody HRP Conjugated |
| Species/Host | Mouse |
| Conjugation | HRP |
| Tested Applications | ICC, IHC |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Clonality | Monoclonal |
| Source | mouse |
| Uniprot ID | P08637 |
| Dilution Range | IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 |
| Expiration Date | 12 months from date of receipt. |