

Product Datasheet

beta I Tubulin(Loading Control) mouse Monoclonal Antibody PE-Cy5 Conjugated (orb1539602)

Description	beta I Tubulin(Loading Control) mouse Monoclonal Antibody PE- Cy5 Conjugated
Species/Host	Mouse
Conjugation	PE/Cy5
Tested Applications	FC, ICC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	mouse
Uniprot ID	Q9H4B7
Dilution Range	Flow-Cyt=1:50-200 ICC=1:50- 200 IF=1:50-200
Expiration Date	12 months from date of receipt.