

Product Datasheet

Myc tag mouse Monoclonal Antibody PE Conjugated (orb1539607)

Description	Myc tag mouse Monoclonal Antibody PE Conjugated
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC, ICC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	mouse
Uniprot ID	P01106
Dilution Range	Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.