

## Product Datasheet

# GFP-Tag mouse Monoclonal Antibody Bio Conjugated (orb1539611)

<b>Description</b>	GFP-Tag mouse Monoclonal Antibody Bio Conjugated
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	WB
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Source</b>	mouse
<b>Uniprot ID</b>	<a href="#">P42212</a>
<b>Dilution Range</b>	WB=1:50-200
<b>Expiration Date</b>	12 months from date of receipt.