

Product Datasheet

phospho-PIK3R1 (Tyr368) Polyclonal Antibody FITC Conjugated (orb1539909)

Description	phospho-PIK3R1 (Tyr368) Polyclonal Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P27986
Dilution Range	IF=1:50-200
Expiration Date	12 months from date of receipt.