

Product Datasheet

MPO Polyclonal Antibody PE-Cy5.5 Conjugated (orb1539951)

Description	MPO Polyclonal Antibody PE-Cy5.5 Conjugated
Species/Host	Rabbit
Conjugation	PE/Cy5.5
Tested Applications	ICC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P05164
Dilution Range	ICC=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.