

Product Datasheet

PRRSV M protein Polyclonal Antibody Bio Conjugated (orb1540025)

Description	PRRSV M protein Polyclonal Antibody Bio Conjugated
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	WB
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q9WJB2
Dilution Range	WB=1:50-200
Expiration Date	12 months from date of receipt.