

Product Datasheet

Phospho-HDAC4 (Ser246) / HDAC5 (Ser259) / HDAC7 (Ser155) Polyclonal Antibody Bio Conjugated (orb1540180)

Description	Phospho-HDAC4 (Ser246) / HDAC5 (Ser259) / HDAC7 (Ser155) Polyclonal Antibody Bio Conjugated
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	IF, IHC, WB
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P56524
Dilution Range	WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200
Expiration Date	12 months from date of receipt.