

Product Datasheet

CD11c Polyclonal Antibody PE-Cy7 Conjugated (orb1540283)

Description	CD11c Polyclonal Antibody PE-Cy7 Conjugated
Species/Host	Rabbit
Conjugation	PE/Cy7
Tested Applications	IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P20702
Dilution Range	IF=1:50-200
Expiration Date	12 months from date of receipt.