

Product Datasheet

Nurr1 Polyclonal Antibody FITC Conjugated (orb1540473)

Description	Nurr1 Polyclonal Antibody FITC Conjugated
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	FC, IF
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	P43354
Dilution Range	Flow-Cyt=2ug/Test IF=1:50-200
Expiration Date	12 months from date of receipt.