

Product Datasheet

PCV2 Cap protein Polyclonal Antibody HRP Conjugated (orb1540566)

Description	PCV2 Cap protein Polyclonal Antibody HRP Conjugated
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	WB
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	O56129
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.