

Product Datasheet

Norovirus Capsid protein VP1 Polyclonal Antibody HRP Conjugated (orb1540603)

Description	Norovirus Capsid protein VP1 Polyclonal Antibody HRP Conjugated
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	ICC, IHC-Fr, IHC-P, WB
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q83884
Dilution Range	WB=1:500-2000 IHC-P=1:50- 200 IHC-F=1:50-200 ICC=1:50- 200
Expiration Date	12 months from date of receipt.