

Product Datasheet

Adenovirus hexon protein Polyclonal Antibody Gold Conjugated (orb1541028)

Description	Adenovirus hexon protein Polyclonal Antibody Gold Conjugated
Species/Host	Rabbit
Conjugation	Gold
Tested Applications	GICA, IEM, IGS
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Recombinant
Source	Rabbit
Uniprot ID	Q9YVE5
Dilution Range	IEM=1:20-200 IGS=1:20-200 GICA=1:20-200
Expiration Date	12 months from date of receipt.