

## Product Datasheet

# PCDH12 Polyclonal Antibody FITC Conjugated (orb1541153)

<b>Description</b>	PCDH12 Polyclonal Antibody FITC Conjugated
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FC, ICC, IF
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Source</b>	Rabbit
<b>Uniprot ID</b>	<a href="#">Q9NPG4</a>
<b>Dilution Range</b>	Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200
<b>Expiration Date</b>	12 months from date of receipt.