

## Product Datasheet

# HLA Class 1 ABC /HLA ABC Polyclonal Antibody PE-Cy5 Conjugated (orb1541211)

<b>Description</b>	HLA Class 1 ABC /HLA ABC Polyclonal Antibody PE-Cy5 Conjugated
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	PE/Cy5
<b>Tested Applications</b>	ICC, IF
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Source</b>	Rabbit
<b>Uniprot ID</b>	<a href="#">P30443</a>
<b>Dilution Range</b>	ICC=1:50-200 IF=1:50-200
<b>Expiration Date</b>	12 months from date of receipt.